





Instructions for VE-11

Handzfree™ Anesthetic Bottle Holder

1. AMI recommends that you have this holder installed by a professional builder or qualified craftsperson. Remove protective covering. Attach the holder to wall using the screws provided. Be sure to position the holder so that the short 'leg' of the L-shaped holder points up and out from the wall.
2. Carefully place anesthetic bottles, necks down, into holder so that anesthetic can be drawn by the practitioner without needing to remove one or more sterile gloves and without need of an assistant.

IMPORTANT: *Walls contain electrical wires, pipes, and other potentially dangerous or easily damaged obstructions. DO NOT* attempt installation unless you have received clearance from your building maintenance personnel, electrician, and other relevant authorities. Also, while the holder is quite light, please check to ensure that the wall can bear its weight and that of 4 full bottles of anesthetic before installing. Advanced Meditech is not responsible or liable for any damage or injury resulting from installation or use of this product.



Thank You For Your Business!

"Technology to Better Health"