Instruments





The cantilever on LiBrand ring clamp is a critical factor in NSV's reduction in complications. As the ring clamp is locked around the vas extracutaneously, its tips spring slightly apart. As a result, they hold firmly without pinching or crushing delicate tissues.

Special Ring Clamps



Pointed Dissector & Others



No-Needle Anesthesia



How Does it Work? Vas

The key to No-Needle Injector is the Madalet. It's a pen-sized, stainless steel injector that uses pneumatic pressure Scrotum to deliver a virtually painless injection of

Nothing except the anesthesic solution penetrates the natients skin. With the push of the release button, anesthesia solution is sprayed through a calibrated orifice. View video of the technige at www.VasectomyStore.com

Marketing Brochures Draw patients to your practice with No-Needle.



anesthesia to vasa

The 2500-V MadaJet includes...

- 1 MadaJet Injection Device
- 1 Full Pack of Madalet Accessories
- 1 No-Needle NSV Manual
- 1 No-Needle NSV DVD Instructional Video
- 1 Sample VED-6NN Marketing Brochures

Cautery Equipment



Geiger TCU Item# 2200 Adjustable Thermal Cautery Unit - 110 Volts

- No unwanted tissue damage Increased cost efficiency
- · No frequency interference

Disposable Thermal Cautery

AM-21 Integrated unit includes handle, vasectomy tip, & AA batteries, Box of 10 units.

Reusable Thermal Cautery

CH-HI Reusable Handle, Burns at 1200°C (2200 °F), nonsterile base unit with 1 nonsterile tip. 1 set replaceable AA batteries.

CT-2121 Sterile Disposable Thermal Cautery Tips for Vasectomy, Individually peel-pouched with dual plastic sheaths to cover non-sterile handle. Box of 10 tips.

Ligating Clips



Item# VE-6-10 Clips

Medium Surgical Clip Closure Systems

4030 Medium Clip Appliers, stainless steel, 6" long, curved tip. For use with 4035 Medium Titanium Surgical Clins only

4035 Medium Surgical Clips, 6 clips per cartridge. Box of 30 cartridges with a total of 180 titanium clins per box. For secure vasal and other occlusion.

VE-6 Medium Clip Appliers, American-made stainless steel, 6" long, curved tip. For use with VE-6-25 and VE-6-10 titanium surgical clips only.

VE-6-10 Medium Surgical Clips, 10 clips per cartridge. Box of 20 cartridges with a total of 200 titanium clips per box. For secure vasal and other occlusion. For use with VE-6 Applier only.

Order at: 800 • 635 • 2452 or visit VasectomyStore.com