


VasectomyStore.com

Your Complete Vasectomy Source Since 1992

Vasal Occlusion

Ligating Clips and Thermal Cautery



Item# 4030DX

The 6-Clip Cartridge Closure System



If you've avoided certain closure systems because you cannot or choose not to resterilize unused clips, this closure system is ideal.

- **Easier to Load**

Widely spaced clips on high contrast cartridge provide superior visibility.

- **Stay-Put Adhesive Cartridge**

Adhesive strip underneath the clip cartridge stabilizes clips for loading.

- **Enhanced Grip**

Clips have positive grip micro-eration design for enhanced grip on vessels.

Hand-held Thermal Cautery

Item# AM-21



Disposable Thermal Cautery

Integrated unit includes handle, vasectomy tip, & AA batteries. Box of 10 units.

Reusable Thermal Cautery & Tips

Burns at approximately 1200°C (2200 °F), non-sterile base unit with 1 non-sterile tip, 2 (AA) replaceable batteries. Additional tips sold separately.

Item# CH-HI



Helpful Materials

The latest studies and articles regarding vasal occlusion are available in the Medical Info section of VasectomyStore.com. Also, find helpful suggestions in the 60 page guide "No-Scalpel Vasectomy - An Illustrated Guide for Surgeons" (item# VED-1) or the DVD "No-Scalpel Vasectomy" (item# VED-2)

